

KERN OBB-A 1544

Centreringsoculair

KERN



Categorie

| | |
|------------------|---------------------|
| Merk | Optics |
| Productcategorie | Optische Instrument |
| Productgroep | Microscoop oculair |

Milieuomstandigheden

| | |
|-------------------------|-------|
| Opslagtemperatuur [Min] | -5 °C |
| Opslagtemperatuur [Max] | 40 °C |

Verpakking & verzending

| | |
|-------------------------------|--------------|
| Afreesbaarheid kracht [d] (N) | 1 d |
| Afmetingen verpakking (B×D×H) | 85×65×45 mm |
| Nettogewicht | 0,2 kg |
| Verzendmethode | Pakketdienst |
| Nettogewicht ca. | 0,20 kg |
| Brutogewicht ca. | 0,40 kg |
| Verzendgewicht | 0,4 kg |

Pictogrammen

STANDARD

