

# KERN OBB-A2406

Eyepiece (Ø 23,2 mm): WF 16 x /Ø 13 mm

# KERN



## Pictograms

### STANDARD



Category	
Brand	Optics
Product category	Optical Instrument
Product group	Microscope eyepiece

Construction	
Anti fungus	✓

Ocular	
Ocular field width	WF
Eye point	Standard
Ocular magnifications	16 x
Ocular visual field	13 mm
Ocular diameter	23,2 mm

Environmental conditions	
Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

Packing & Shipping	
Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	85×65×45 mm
Net weight	0,1 kg
Shipping method	Parcel service
Net weight approx.	0,10 kg
Gross weight approx.	0,40 kg
Shipping weight	0,4 kg