

KERN OBB-A2424

Fixed Stage for inverted microscopes

KERN



Category

Brand	Optics
Product category	Optical Instrument
Product group	Microscope stage

Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	280×260×40 mm
Net weight	2,5 kg
Shipping method	Parcel service
Net weight approx.	2,6 kg
Gross weight approx.	3,0 kg
Shipping weight	3 kg

Pictograms