

KERN ODC-A2403

Calibration specimen slide with micrometer (0,01 mm)

KERN



Category

Brand	Optics
Product category	Optical Instrument
Product group	Calibration Specimen Slide

Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	80×50×40 mm
Net weight	0,1 kg
Shipping method	Parcel service
Net weight approx.	0,10 kg
Gross weight approx.	0,20 kg
Shipping weight	0,2 kg

Pictograms

STANDARD

