

KERN ORA 5PB

KERN

Refractive index measurement for laboratories and the industry



Construction

Dimension (W×D×H)	160×40×40 mm
Refractometer type	Handheld refractometer
Scale - Display type	Analogue
Material prism	ZLaF3 lanthanum flint glass nd=1.85544

Functions

IP protection - complete device -

Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

Packing & Shipping

Dimensions transport container (W×D×H)	205×75×55 mm
Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	215×94×65 mm
Net weight	0,30372 kg
Shipping method	Parcel service
Net weight approx.	0,35 kg
Gross weight approx.	0,40 kg
Shipping weight	0,356 kg

Services (optional)

Article number for factory calibration	961-290
--	---------

Pictograms

STANDARD



Category

Brand	Optics
Product category	Refractometer
Product group	Analogue refractometer
Product family	ORA-P

Measuring System

Measuring method	Light refraction
Area of application	Urine
Appropriate calibration liquid	0 % (distilled water)
Measurement temperature [Min] (°C)	10 °C
Measurement temperature [Max] (°C)	20 °C
Scale protein serum	✓
Readout protein serum	0,1 g/100 ml
Measuring range protein serum [Min]	2 g/100 ml
Measuring range protein serum [Max]	14 g/100 ml
Scale urine (s. g. dog)	✓
Readout urine (s. g. dog)	0,001
Measuring range urine (s. g. dog) [Min]	1,000
Measuring range urine (s. g. dog) [Max]	1,060
Scale urine (s. g. cat)	✓
Readout urine (s. g. cat)	0,001
Measuring range urine (s. g. cat) [Min]	1,000
Measuring range urine (s. g. cat) [Max]	1,060