

KERN OZB-A4207

Anti-fungal 0,5 x

KERN



Category

Brand	Optics
Product category	Optical Instrument
Product group	Microscope objective lens

Construction

Anti fungus	✓
-------------	---

Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	55×55×25 mm
Net weight	0,7 kg
Shipping method	Parcel service
Net weight approx.	0,70 kg
Gross weight approx.	1,00 kg
Shipping weight	1 kg

Objective

Objective magnification	0,5 x
Objective working distance	250 mm

Pictograms

STANDARD

