

KERN OZB-A4532



Double gooseneck 490 mm; in conjunction with OZB-A4512



Category

Brand	Optics
Product category	Optical Instrument
Product group	Microscope illumination

Construction

Beleuchtungsart	Halogen
Gooseneck length	490 mm

Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	280×185×50 mm
Net weight	0,7 kg
Shipping method	Parcel service
Net weight approx.	0,70 kg
Gross weight approx.	1,00 kg
Shipping weight	1 kg

Pictograms

STANDARD

