

# KERN OZB-A4601

Dark field kit for stereo microscopes

**KERN**



## Category

Brand	Optics
Product category	Optical Instrument
Product group	Dark field inserts

## Construction

Material	metal
----------	-------

## Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	125×125×45 mm
Net weight	0,7 kg
Shipping method	Parcel service
Net weight approx.	0,70 kg
Gross weight approx.	1,00 kg
Shipping weight	1 kg

## Pictograms

STANDARD

