

# KERN OZB-A5203

**KERN**

Individuality, variety and flexible working through our modular construction system ► PREMIUM-Universal stands



## Category

Brand	Optics
Product category	Microscope
Product group	Stereo Microscope- Stand
Product family	OZB-UP

## Construction

Dimension (W×D×H)	254×545×515 mm
-------------------	----------------

## Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	680×330×260 mm
Net weight	21,1772 kg
Shipping method	Parcel service
Net weight approx.	22 kg
Gross weight approx.	24 kg
Shipping weight	22,5 kg

## Pictograms

STANDARD

