

KERN OZB-A56 14

KERN

Gesoldeerde bescherm lens voor stereomicroscopen



Categorie

Merk	Optics
Productcategorie	Optische Instrument
Productgroep	Lens

Milieuomstandigheden

Opslagtemperatuur [Min]	-5 °C
Opslagtemperatuur [Max]	40 °C

Verpakking & verzending

Afreesbaarheid kracht [d] (N)	1 d
Afmetingen verpakking (B×D×H)	55×55×25 mm
Nettogewicht	0,5 kg
Verzendmethode	Pakketdienst
Nettogewicht ca.	0,50 kg
Brutogewicht ca.	0,70 kg
Verzendgewicht	0,7 kg

Pictogrammen

STANDARD

