

# KERN OZB-A7101EU

**KERN**

Professional illumination units guarantee outstanding, uniform and strong illumination



Gross weight approx. 1,00 kg  
Shipping weight 1 kg

## Pictograms

### STANDARD



## Category

Brand	Optics
Product category	Microscope
Product group	Microscope illumination
Product family	OZB-IR

## Approval

CE mark	✓
---------	---

## Construction

Dimension (W×D×H)	99×90×33 mm
Beleuchtungsart	LED
Illuminance	4,5 W
Colour temperature	6500-7000 K
Illumination dimmable	✓
Ring illumination diameter	6,2 cm

## Power Supply

Supplied power supply	Power supply unit
Plug-in power supply / adapter for countries - optional	EURO

## Environmental conditions

Storage temperature [Min]	0 °C
Storage temperature [Max]	40 °C

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	250×150×80 mm
Net weight	0,9375 kg
Shipping method	Parcel service
Net weight approx.	0,95 kg