

KERN 307-01

Individual weight, OIML Class E1, 1 g, Knob, Stainless steel polished



Services (optional)

Article number for DAkkS calibration	963-231
Article number for guaranteed plus-deviation	921-231

Pictograms

OPTION

DAkkS
+3 DAYS

Category

Brand	KERN
Product category	Individual weight
Product group	OIML E1
Product family	307

Measuring System

Test weight - nominal value	1 g
Test weight - volume (cm ³)	0,125
Test weight - volume uncertainty (cm ³)	0,0016
Permissible tolerance according to OIML (+/-)	0,01 mg
Test weight - Density [Min] (OIML)	6900
Test weight - Density [Max] (OIML)	9600
Weight OIML class	E1

Construction

Shape weight	Knob
Dimensions (Ø×H)	6×6 mm
Material	Stainless steel polished

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	6×6×6 mm
Net weight	0,001 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg