

# KERN 307-05

**KERN**

Individual weight, OIML Class E1, 20 g, Knob, Stainless steel polished



## Services (optional)

Article number for DAkkS calibration	963-235
Article number for guaranteed plus-deviation	921-235

## Pictograms

### OPTION

DAkkS  
+3 DAYS

## Category

Brand	KERN
Product category	Individual weight
Product group	OIML E1
Product family	307

## Measuring System

Test weight - nominal value	20 g
Test weight - volume (cm <sup>3</sup> )	2,5
Test weight - volume uncertainty (cm <sup>3</sup> )	0,0313
Permissible tolerance according to OIML (+/-)	0,025 mg
Test weight - Density [Min] (OIML)	7840
Test weight - Density [Max] (OIML)	8170
Weight OIML class	E1

## Construction

Shape weight	Knob
Dimensions (Ø×H)	13×22 mm
Material	Stainless steel polished

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	13×13×22 mm
Net weight	0,02 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg