

KERN 307-06

KERN

Individual weight, OIML Class E1, 50 g, Knob, Stainless steel polished



Services (optional)

Article number for DAkkS calibration	963-236
Article number for guaranteed plus-deviation	921-236

Pictograms

OPTION

DAkkS
+3 DAYS

Category

Brand	KERN
Product category	Individual weight
Product group	OIML E1
Product family	307

Measuring System

Test weight - nominal value	50 g
Test weight - volume (cm ³)	6,25
Test weight - volume uncertainty (cm ³)	0,0781
Permissible tolerance according to OIML (+/-)	0,03 mg
Test weight - Density [Min] (OIML)	7920
Test weight - Density [Max] (OIML)	8080
Weight OIML class	E1

Construction

Shape weight	Knob
Dimensions (Ø×H)	18×28 mm
Material	Stainless steel polished

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	18×18×28 mm
Net weight	0,05 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg