

# KERN 307-07

**KERN**

Individual weight, OIML Class E1, 100 g, Knob, Stainless steel polished



## Services (optional)

|  |         |
|--|---------|
| Article number for DAkkS calibration         | 963-237 |
| Article number for guaranteed plus-deviation | 921-237 |

## Pictograms

## OPTION

DAkkS  
+3 DAYS

## Category

|                  |                   |
|------------------|-------------------|
| Brand            | KERN              |
| Product category | Individual weight |
| Product group    | OIML E1           |
| Product family   | 307               |

## Measuring System

|   |         |
|---|---------|
| Test weight - nominal value                         | 100 g   |
| Test weight - volume (cm <sup>3</sup> )             | 12,5    |
| Test weight - volume uncertainty (cm <sup>3</sup> ) | 0,1563  |
| Permissible tolerance according to OIML (+/-)       | 0,05 mg |
| Test weight - Density [Min] (OIML)                  | 7934    |
| Test weight - Density [Max] (OIML)                  | 8067    |
| Weight OIML class                                   | E1      |

## Construction

|                  |                          |
|------------------|--------------------------|
| Shape weight     | Knob                     |
| Dimensions (Ø×H) | 22×38 mm                 |
| Material         | Stainless steel polished |

## Packing & Shipping

|                              |                |
|------------------------------|----------------|
| Readability force [d] (N)    | 1 d            |
| Dimensions packaging (W×D×H) | 22×22×38 mm    |
| Net weight                   | 0,1 kg         |
| Shipping method              | Parcel service |
| Net weight approx.           | 0,10 kg        |
| Gross weight approx.         | 0,10 kg        |
| Shipping weight              | 0,1 kg         |