

KERN 317-01

Individual weight, OIML Class E2, 1 g, Knob, stainless steel polished (OIML)



Services (optional)

Article number for conformity assessment	952-331
Article number for DAkkS calibration	962-331
Article number for guaranteed plus-deviation	921-331

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML E2
Product family	317

Measuring System

Test weight - nominal value	1 g
Test weight - volume (cm ³)	0,125
Test weight - volume uncertainty (cm ³)	0,0016
Permissible tolerance according to OIML (+/-)	0,03 mg
Test weight - Density [Min] (OIML)	5300
Test weight - Density [Max] (OIML)	16000
Weight OIML class	E2

Construction

Shape weight	Knob
Dimensions (Ø×H)	6×6 mm
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	6×6×6 mm
Net weight	0,001 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg