

# KERN 317-02

Individual weight, OIML Class E2, 2 g, Knob, stainless steel polished (OIML)



## Services (optional)

Article number for conformity assessment	952-332
Article number for DAkkS calibration	962-332
Article number for guaranteed plus-deviation	921-332

## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML E2
Product family	317

## Measuring System

Test weight - nominal value	2 g
Test weight - volume (cm <sup>3</sup> )	0,25
Test weight - volume uncertainty (cm <sup>3</sup> )	0,0031
Permissible tolerance according to OIML (+/-)	0,04 mg
Test weight - Density [Min] (OIML)	6000
Test weight - Density [Max] (OIML)	12000
Weight OIML class	E2

## Construction

Shape weight	Knob
Dimensions (Ø×H)	6×10 mm
Material	stainless steel polished (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	6×6×10 mm
Net weight	0,002 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg