

KERN 317-06

Individual weight, OIML Class E2, 50 g, Knob, stainless steel polished (OIML)



Services (optional)

Article number for conformity assessment	952-336
Article number for DAkkS calibration	962-336
Article number for guaranteed plus-deviation	921-336

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML E2
Product family	317

Measuring System

Test weight - nominal value	50 g
Test weight - volume (cm ³)	6,25
Test weight - volume uncertainty (cm ³)	0,0781
Permissible tolerance according to OIML (+/-)	0,1 mg
Test weight - Density [Min] (OIML)	7740
Test weight - Density [Max] (OIML)	8890
Weight OIML class	E2

Construction

Shape weight	Knob
Dimensions (Ø×H)	18×28 mm
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	18×18×28 mm
Net weight	0,05 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg