

KERN 317-07

KERN

Individual weight, OIML Class E2, 100 g, Knob, stainless steel polished (OIML)



Services (optional)

Article number for conformity assessment	952-336
Article number for DAkkS calibration	962-337
Article number for guaranteed plus-deviation	921-337

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML E2
Product family	317

Measuring System

Test weight - nominal value	100 g
Test weight - volume (cm ³)	12,5
Test weight - volume uncertainty (cm ³)	0,1563
Permissible tolerance according to OIML (+/-)	0,16 mg
Test weight - Density [Min] (OIML)	7810
Test weight - Density [Max] (OIML)	8210
Weight OIML class	E2

Construction

Shape weight	Knob
Dimensions (Ø×H)	22×38 mm
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	22×22×38 mm
Net weight	0,1 kg
Shipping method	Parcel service
Net weight approx.	0,10 kg
Gross weight approx.	0,10 kg
Shipping weight	0,1 kg