

# KERN 328-01

Milligram weights, F1, 1 mg, Platelet, stainless steel polished (OIML)



## Services (optional)

Article number for conformity assessment	952-451
Article number for DAkkS calibration	962-451
Article number for guaranteed plus-deviation	921-451

## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	328

## Measuring System

Test weight - nominal value	1 mg
Test weight - volume (cm <sup>3</sup> )	0,0004
Test weight - volume uncertainty (cm <sup>3</sup> )	0,0001
Permissible tolerance according to OIML (+/-)	0,02 mg
Weight OIML class	F1

## Construction

Shape weight	Platelet
Dimensions weight (W×D×H)	3,5×3,5×3 mm
Material	stainless steel polished (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	3,5×3,5×3 mm
Net weight	0,01 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg