

# KERN 346-76

Individual weight, OIML Class M1, 5 kg, Block Eco, cast iron lacquered



## Services (optional)

Article number for conformity assessment	952-643
Article number for DAkkS calibration	962-643

## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	346

## Measuring System

Test weight - nominal value	5 kg
Test weight - volume (cm <sup>3</sup> )	704,2254
Test weight - volume uncertainty (cm <sup>3</sup> )	59,512
Permissible tolerance according to OIML (+/-)	250 mg
Test weight - Density [Min] (OIML)	3000
Weight OIML class	M1

## Construction

Shape weight	Block Eco
Dimensions weight (W×D×H)	150×79×89 mm
Material	cast iron lacquered

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	170×95×100 mm
Net weight	5 kg
Shipping method	Parcel service
Net weight approx.	5,0 kg
Gross weight approx.	5,0 kg
Shipping weight	5 kg