

KERN 346-77

Individual weight, OIML Class M1, 10 kg, Block Eco, cast iron lacquered



Services (optional)

Article number for conformity assessment	952-644
Article number for DAkkS calibration	962-644

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	346

Measuring System

Test weight - nominal value	10 kg
Test weight - volume (cm ³)	1408,4507
Test weight - volume uncertainty (cm ³)	119,024
Permissible tolerance according to OIML (+/-)	500 mg
Test weight - Density [Min] (OIML)	3000
Weight OIML class	M1

Construction

Shape weight	Block Eco
Dimensions weight (W×D×H)	190×95×109 mm
Material	cast iron lacquered

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	210×110×125 mm
Net weight	10 kg
Shipping method	Parcel service
Net weight approx.	10 kg
Gross weight approx.	10 kg
Shipping weight	10 kg