


# KERN 346-89

**KERN**

Individual weight, OIML Class M1, 50 kg, Block, cast iron lacquered



## Services (optional)

Article number for conformity assessment	952-646
Article number for DAkkS calibration	962-646

## Pictograms

## STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	346

## Measuring System

Test weight - nominal value	50 kg
Test weight - volume (cm <sup>3</sup> )	7042,2535
Test weight - volume uncertainty (cm <sup>3</sup> )	595,12
Permissible tolerance according to OIML (+/-)	2,5 g
Test weight - Density [Min] (OIML)	3000
Weight OIML class	M1

## Construction

Shape weight	Block
Dimensions weight (W×D×H)	310×155×192 mm
Material	cast iron lacquered

## Packing & Shipping

Readability force [d] (N)	2 d
Dimensions packaging (W×D×H)	330×180×230 mm
Net weight	50 kg
Shipping method	Freight
Net weight approx.	50 kg
Gross weight approx.	50 kg
Shipping weight	50 kg