

KERN 347-015

KERN

Individual weight, OIML Class M1, 1 g, Slot weight (Type A), stainless steelfine turned (OIML)



Services (optional)

Article number for DAkkS calibration 962-631

Pictograms

STANDARD



Category	
Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	347

Measuring System	
Test weight - nominal value	1 g
Test weight - volume (cm ³)	0,119
Test weight - volume uncertainty (cm ³)	0,0024
Permissible tolerance according to OIML (+/-)	1 mg
Test weight - Density [Min] (OIML)	2600
Weight OIML class	M1

Construction	
Shape weight	Slot weight (Type A)
Dimensions (Ø×H)	13×2 mm
Material	stainless steelfine turned (OIML)

Packing & Shipping	
Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	13×13×2 mm
Net weight	0,001 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg