

# KERN 347-015-N00

**KERN**

Individual weight (Newton), 1 cN, Slot weight (Type A), stainless steelfine turned



## Category

Brand	KERN
Product category	Individual weight
Product group	Newton weight
Product family	347-N

## Measuring System

Test weight - nominal value (N)	1 cN
---------------------------------	------

## Construction

Shape weight	Slot weight (Type A)
Material	stainless steelfine turned
Dimensions (Ø×H)	13×2 mm

## Packing & Shipping

Shipping method	Parcel service
Delivery time	1 d
Dimensions packaging (W×D×H)	13×13×2 mm
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

## Pictograms

STANDARD



1 DAY