

KERN 347-035-N00

KERN

Individual weight (Newton), 5 cN, Slot weight (Type A), stainless steelfine turned



Category

Brand	KERN
Product category	Individual weight
Product group	Newton weight
Product family	347-N

Measuring System

Test weight - nominal value (N)	5 cN
---------------------------------	------

Construction

Shape weight	Slot weight (Type A)
Material	stainless steelfine turned
Dimensions (Ø×H)	22×5 mm

Packing & Shipping

Shipping method	Parcel service
Delivery time	1 d
Dimensions packaging (W×D×H)	18×18×4 mm
Net weight approx.	0,05 kg
Gross weight approx.	0,10 kg
Shipping weight	0,1 kg

Pictograms

STANDARD

