

# KERN 347-036

Individual weight, OIML Class M1, 5 g, Hook weight, stainless steelfine turned (OIML)



## Services (optional)

Article number for DAkkS calibration	962-633
--------------------------------------	---------

## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	347

## Measuring System

Test weight - nominal value	5 g
Test weight - volume (cm <sup>3</sup> )	0,5952
Test weight - volume uncertainty (cm <sup>3</sup> )	0,012
Permissible tolerance according to OIML (+/-)	1,6 mg
Test weight - Density [Min] (OIML)	2600
Weight OIML class	M1

## Construction

Shape weight	Hook weight
Dimensions (Ø×H)	8×22 mm
Material	stainless steelfine turned (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	8×8×22 mm
Net weight	0,005 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg