

KERN 347-045

KERN

Individual weight, OIML Class M1, 10 g, Slot weight (Type A), stainless steelfine turned (OIML)



Services (optional)

Article number for DAkkS calibration 962-634

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	347

Measuring System

Test weight - nominal value	10 g
Test weight - volume (cm ³)	1,1905
Test weight - volume uncertainty (cm ³)	0,0241
Permissible tolerance according to OIML (+/-)	2 mg
Test weight - Density [Min] (OIML)	2600
Weight OIML class	M1

Construction

Shape weight	Slot weight (Type A)
Dimensions (Ø×H)	22×5 mm
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	22×22×5 mm
Net weight	0,01 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg