

# KERN 347-047

**KERN**

Individual weight, OIML Class M1, 10 g, Slot weight (Type B), stainless steelfine turned (OIML)



## Pictograms

### STANDARD



### Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	347

### Measuring System

Test weight - nominal value	10 g
Test weight - volume (cm <sup>3</sup> )	1,1905
Test weight - volume uncertainty (cm <sup>3</sup> )	0,0241
Permissible tolerance according to OIML (+/-)	2 mg
Test weight - Density [Min] (OIML)	2600
Weight OIML class	M1

### Construction

Shape weight	Slot weight (Type B)
Dimensions (Ø×H)	28×5 mm
Material	stainless steelfine turned (OIML)

### Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	28×28×5 mm
Net weight	0,01 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

### Services (optional)

Article number for DAkkS calibration	962-634
--------------------------------------	---------