

KERN 347-055-N00

KERN

Individual weight (Newton), 20 cN, Slot weight (Type A), stainless steelfine turned



Category

Brand	KERN
Product category	Individual weight
Product group	Newton weight
Product family	347-N

Measuring System

Test weight - nominal value (N)	20 cN
---------------------------------	-------

Construction

Shape weight	Slot weight (Type A)
Material	stainless steelfine turned
Dimensions (Ø×H)	27×6 mm

Packing & Shipping

Shipping method	Parcel service
Delivery time	1 d
Dimensions packaging (W×D×H)	27×27×6 mm
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

Pictograms

STANDARD

