

KERN 347-056-N00

KERN

Individual weight (Newton), 20 cN, Hook weight, stainless steelfine turned (OIML)



Category

Brand	KERN
Product category	Individual weight
Product group	Newton weight
Product family	347-N

Measuring System

Test weight - nominal value (N)	20 cN
---------------------------------	-------

Construction

Shape weight	Hook weight
Material	stainless steelfine turned (OIML)
Dimensions (Ø×H)	13×30 mm

Packing & Shipping

Shipping method	Parcel service
Delivery time	1 d
Dimensions packaging (W×D×H)	13×13×30 mm
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

Pictograms

STANDARD

