

KERN 347-057

KERN

Individual weight, OIML Class M1, 20 g, Slot weight (Type B), stainless steelfine turned (OIML)



Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	347

Measuring System

Test weight - nominal value	20 g
Test weight - volume (cm ³)	2,381
Test weight - volume uncertainty (cm ³)	0,0482
Permissible tolerance according to OIML (+/-)	2,5 mg
Test weight - Density [Min] (OIML)	2600
Weight OIML class	M1

Construction

Shape weight	Slot weight (Type B)
Dimensions (Ø×H)	28×3,8 mm
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	28×28×4 mm
Net weight	0,02 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

Services (optional)

Article number for DAkkS calibration	962-635
--------------------------------------	---------