

# KERN 347-057-N00

**KERN**

Individual weight (Newton), 20 cN, Slot weight (Type B), stainless steelfine turned



## Category

Brand	KERN
Product category	Individual weight
Product group	Newton weight
Product family	347-N

## Measuring System

Test weight - nominal value (N)	20 cN
---------------------------------	-------

## Construction

Shape weight	Slot weight (Type B)
Material	stainless steelfine turned
Dimensions (Ø×H)	38×3,95 mm

## Packing & Shipping

Shipping method	Parcel service
Delivery time	1 d
Dimensions packaging (W×D×H)	38×38×4 mm
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

## Pictograms

STANDARD

