

KERN 347-075-100

KERN

Supporting bar, OIML Class M1, 100 g, stainless steelfine turned (OIML)



Services (optional)

Article number for DAkkS calibration	962-637
--------------------------------------	---------

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	Supporting Bar
Product family	347

Measuring System

Test weight - nominal value	100 g
Test weight - volume (cm ³)	11,9048
Test weight - volume uncertainty (cm ³)	0,2409
Permissible tolerance according to OIML (+/-)	5 mg
Test weight - Density [Min] (OIML)	2600
Weight OIML class	M1

Construction

Shape weight	Supporting bar (Type A)
Max. total weight	2 kg
Max. single weight	1 kg
Dimensions (Ø×H)	26×240 mm
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	26×26×240 mm
Net weight	0,1 kg
Shipping method	Parcel service
Net weight approx.	0,10 kg
Gross weight approx.	0,10 kg
Shipping weight	0,1 kg