

# KERN 347-076

# KERN

Individual weight, OIML Class M1, 100 g, Hook weight, stainless steelfine turned (OIML)



## Services (optional)

Article number for DAkkS calibration	962-637
--------------------------------------	---------

## Pictograms

## STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	347

## Measuring System

Test weight - nominal value	100 g
Test weight - volume (cm <sup>3</sup> )	11,9048
Test weight - volume uncertainty (cm <sup>3</sup> )	0,2409
Permissible tolerance according to OIML (+/-)	5 mg
Test weight - Density [Min] (OIML)	2600
Weight OIML class	M1

## Construction

Shape weight	Hook weight
Dimensions (Ø×H)	22×54 mm
Material	stainless steelfine turned (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	22×22×54 mm
Net weight	0,1 kg
Shipping method	Parcel service
Net weight approx.	0,10 kg
Gross weight approx.	0,10 kg
Shipping weight	0,1 kg