

# KERN 347-086

# KERN

Individual weight, OIML Class M1, 200 g, Hook weight, stainless steelfine turned (OIML)



## Services (optional)

Article number for DAkkS calibration 962-638

## Pictograms

## STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	347

## Measuring System

Test weight - nominal value	200 g
Test weight - volume (cm <sup>3</sup> )	23,8095
Test weight - volume uncertainty (cm <sup>3</sup> )	0,4819
Permissible tolerance according to OIML (+/-)	10 mg
Test weight - Density [Min] (OIML)	2600
Weight OIML class	M1

## Construction

Shape weight	Hook weight
Dimensions (Ø×H)	28×61 mm
Material	stainless steelfine turned (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	28×28×61 mm
Net weight	0,2 kg
Shipping method	Parcel service
Net weight approx.	0,20 kg
Gross weight approx.	0,20 kg
Shipping weight	0,2 kg