


KERN 347-087-N00

KERN

Individual weight (Newton), 2 N, Slot weight (Type B), stainless steelfine turned



Category

Brand	KERN
Product category	Individual weight
Product group	Newton weight
Product family	347-N

Measuring System

Test weight - nominal value (N)	2 N
---------------------------------	-----

Construction

Shape weight	Slot weight (Type B)
Material	stainless steelfine turned
Dimensions (Ø×H)	60×14,1 mm

Packing & Shipping

Shipping method	Parcel service
Delivery time	1 d
Dimensions packaging (W×D×H)	60×60×14,1 mm
Net weight approx.	0,20 kg
Gross weight approx.	0,20 kg
Shipping weight	0,2 kg

Pictograms

STANDARD

