

KERN 347-095-100

KERN

Supporting bar, OIML Class M1, 500 g, stainless steelfine turned (OIML)



Services (optional)

Article number for DAkkS calibration 962-639

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	Supporting Bar
Product family	347

Measuring System

Test weight - nominal value	500 g
Test weight - volume (cm ³)	59,5238
Test weight - volume uncertainty (cm ³)	1,2046
Permissible tolerance according to OIML (+/-)	25 mg
Test weight - Density [Min] (OIML)	2600
Weight OIML class	M1

Construction

Shape weight	Supporting bar (Type A)
Max. total weight	20 kg
Max. single weight	10 kg
Dimensions (Ø×H)	40×663 mm
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	40×40×663 mm
Net weight	0,5 kg
Shipping method	Parcel service
Net weight approx.	0,50 kg
Gross weight approx.	0,50 kg
Shipping weight	0,5 kg