

KERN 347-097-N00

KERN

Individual weight (Newton), 5 N, Slot weight (Type B), stainless steelfine turned



Category

| | |
|------------------|-------------------|
| Brand | KERN |
| Product category | Individual weight |
| Product group | Newton weight |
| Product family | 347-N |

Measuring System

| | |
|---------------------------------|-----|
| Test weight - nominal value (N) | 5 N |
|---------------------------------|-----|

Construction

| | |
|------------------|----------------------------|
| Shape weight | Slot weight (Type B) |
| Material | stainless steelfine turned |
| Dimensions (Ø×H) | 60×31,9 mm |

Packing & Shipping

| | |
|------------------------------|----------------|
| Shipping method | Parcel service |
| Delivery time | 1 d |
| Dimensions packaging (W×D×H) | 60×60×32 mm |
| Net weight approx. | 0,50 kg |
| Gross weight approx. | 0,50 kg |
| Shipping weight | 0,5 kg |

Pictograms

STANDARD

