

# KERN 347-115-100

Supporting bar, OIML Class M1, 1 kg, stainless steelfine turned (OIML)



## Services (optional)

Article number for DAkkS calibration	962-641
--------------------------------------	---------

## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	Supporting Bar
Product family	347

## Measuring System

Test weight - nominal value	1 kg
Test weight - volume (cm <sup>3</sup> )	119,0476
Test weight - volume uncertainty (cm <sup>3</sup> )	2,4093
Permissible tolerance according to OIML (+/-)	50 mg
Test weight - Density [Min] (OIML)	2600
Weight OIML class	M1

## Construction

Shape weight	Supporting bar (Type A)
Max. total weight	40 kg
Max. single weight	10 kg
Dimensions (Ø×H)	60×991 mm
Material	stainless steelfine turned (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	60×60×991 mm
Net weight	1 kg
Shipping method	Parcel service
Net weight approx.	1,00 kg
Gross weight approx.	1,00 kg
Shipping weight	1 kg