

KERN 347-117-100-N00

KERN

Supporting bar (Newton), 10 N, stainless steelfine turned



Category

Brand	KERN
Product category	Individual weight
Product group	Supporting Bar (Newton)
Product family	347-N

Measuring System

Test weight - nominal value (N)	10 N
---------------------------------	------

Construction

Shape weight	Supporting bar (Type B)
Material	stainless steelfine turned
Max. total weight (Newton)	460 N
Max. single weight (Newton)	100 N
Dimensions (Ø×H)	55×384 mm

Packing & Shipping

Shipping method	Parcel service
Delivery time	1 d
Dimensions packaging (W×D×H)	55×55×384 mm
Net weight approx.	1,00 kg
Gross weight approx.	1,00 kg
Shipping weight	1 kg

Pictograms

STANDARD

