

KERN 347-117-N00

KERN

Individual weight (Newton), 10 N, Slot weight (Type B), stainless steelfine turned



Category

Brand	KERN
Product category	Individual weight
Product group	Newton weight
Product family	347-N

Measuring System

Test weight - nominal value (N)	10 N
---------------------------------	------

Construction

Shape weight	Slot weight (Type B)
Material	stainless steelfine turned
Dimensions (Ø×H)	118×17,1 mm

Packing & Shipping

Shipping method	Parcel service
Delivery time	1 d
Dimensions packaging (W×D×H)	118×118×17,1 mm
Net weight approx.	1,00 kg
Gross weight approx.	1,00 kg
Shipping weight	1 kg

Pictograms

STANDARD

