

KERN 347-137

KERN

Individual weight, OIML Class M1, 5 kg, Slot weight (Type B), stainless steelfine turned (OIML)



Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	347

Measuring System

Test weight - nominal value	5 kg
Test weight - volume (cm ³)	595,2381
Test weight - volume uncertainty (cm ³)	12,0465
Permissible tolerance according to OIML (+/-)	250 mg
Test weight - Density [Min] (OIML)	2600
Weight OIML class	M1

Construction

Shape weight	Slot weight (Type B)
Dimensions (Ø×H)	178×33,6 mm
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	178×178×34 mm
Net weight	5 kg
Shipping method	Parcel service
Net weight approx.	5,0 kg
Gross weight approx.	5,0 kg
Shipping weight	5 kg

Services (optional)

Article number for DAkkS calibration	962-643
--------------------------------------	---------