

KERN 347-147-100-N00

KERN

Supporting bar (Newton), 100 N, stainless steelfine turned



Category

Brand	KERN
Product category	Individual weight
Product group	Supporting Bar (Newton)
Product family	347-N

Measuring System

Test weight - nominal value (N)	100 N
---------------------------------	-------

Construction

Shape weight	Supporting bar (Type B)
Material	stainless steelfine turned
Max. total weight (Newton)	1,26 kN
Max. single weight (Newton)	100 N
Dimensions (Ø×H)	178×901,9 mm

Packing & Shipping

Shipping method	Parcel service
Delivery time	1 d
Dimensions packaging (W×D×H)	178×178×901,9 mm
Net weight approx.	10 kg
Gross weight approx.	10 kg
Shipping weight	10 kg

Pictograms

STANDARD

