


# KERN 347-147

# KERN

Individual weight, OIML Class M1, 10 kg, Slot weight (Type B), stainless steelfine turned (OIML)



## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	347

## Measuring System

Test weight - nominal value	10 kg
Test weight - volume (cm <sup>3</sup> )	1190,4762
Test weight - volume uncertainty (cm <sup>3</sup> )	24,093
Permissible tolerance according to OIML (+/-)	500 mg
Test weight - Density [Min] (OIML)	2600
Weight OIML class	M1

## Construction

Shape weight	Slot weight (Type B)
Dimensions (Ø×H)	178×61,1 mm
Material	stainless steelfine turned (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	178×178×62 mm
Net weight	10 kg
Shipping method	Parcel service
Net weight approx.	10 kg
Gross weight approx.	10 kg
Shipping weight	10 kg

## Services (optional)

Article number for DAkkS calibration	962-644
--------------------------------------	---------