


KERN 347-147-N00

KERN

Individual weight (Newton), 100 N, Slot weight (Type B), stainless steelfine turned



Category

Brand	KERN
Product category	Individual weight
Product group	Newton weight
Product family	347-N

Measuring System

Test weight - nominal value (N)	100 N
---------------------------------	-------

Construction

Shape weight	Slot weight (Type B)
Material	stainless steelfine turned
Dimensions (Ø×H)	178×61,1 mm

Packing & Shipping

Shipping method	Parcel service
Delivery time	1 d
Dimensions packaging (W×D×H)	178×178×62 mm
Net weight approx.	10 kg
Gross weight approx.	10 kg
Shipping weight	10 kg

Pictograms

STANDARD

