

# KERN 348-05

# KERN

Milligram weights, M1, 20 mg, Platelet, stainless steel polished (OIML)



## Pictograms

### STANDARD



### Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	348

### Measuring System

Test weight - nominal value	20 mg
Permissible tolerance according to OIML (+/-)	0,3 mg
Weight OIML class	M1

### Construction

Shape weight	Platelet
Dimensions weight (W×D×H)	5,5×5,5×2 mm
Material	stainless steel polished (OIML)

### Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	5,5×5,5×2 mm
Net weight	0,01 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

### Services (optional)

Article number for conformity assessment	952-655
Article number for DAkkS calibration	962-655