

KERN 348-07

KERN

Milligram weights, M1, 100 mg, Platelet, stainless steel polished (OIML)



Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	348

Measuring System

Test weight - nominal value	100 mg
Permissible tolerance according to OIML (+/-)	0,5 mg
Weight OIML class	M1

Construction

Shape weight	Platelet
Dimensions weight (W×D×H)	8,5×8,5×4 mm
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	8,5×8,5×4 mm
Net weight	0,0001 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

Services (optional)

Article number for conformity assessment	952-657
Article number for DAkkS calibration	962-657