

KERN 356-86

KERN

Individual weight, OIML Class M2, 5 kg, Block, cast iron lacquered



Services (optional)

Article number for conformity assessment	952-643
Article number for DAkkS calibration	962-643

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M2
Product family	356

Measuring System

Test weight - nominal value	5 kg
Test weight - volume (cm ³)	704,2254
Test weight - volume uncertainty (cm ³)	59,512
Permissible tolerance according to OIML (+/-)	800 mg
Test weight - Density [Min] (OIML)	2300
Weight OIML class	M2

Construction

Shape weight	Block
Dimensions weight (W×D×H)	150×75×84 mm
Material	cast iron lacquered

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	150×75×84 mm
Net weight	5 kg
Shipping method	Parcel service
Net weight approx.	5,0 kg
Gross weight approx.	5,0 kg
Shipping weight	5 kg