

KERN 357-03

Individual weight, OIML Class M2, 5 g, Knob, stainless steelfine turned (OIML)



Services (optional)

Article number for DAkkS calibration	962-633
--------------------------------------	---------

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M2
Product family	357

Measuring System

Test weight - nominal value	5 g
Test weight - volume (cm ³)	0,5952
Test weight - volume uncertainty (cm ³)	0,012
Permissible tolerance according to OIML (+/-)	5 mg
Weight OIML class	M2

Construction

Shape weight	Knob
Dimensions (Ø×H)	8×14 mm
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	24×24×51 mm
Net weight	0,005 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg